

ORBITAL SHAKER

SK-O SERIES

SK-O series orbital shaker is easy and convenient to operate, with stable and low-energy consumption, suitable for mixing different samples. Different clamping methods can easily handle the mixing of samples in commonly used containers such as culture dishes, culture bottles, beakers, flasks, conical flasks, etc.



PRODUCT FEATURES

- High speed control precision, and two operating modes can be selected: timed operation or continuous operation. The positive timed operation time can be displayed in continuous operation (limited to T-shaped timed models);
- DC brushless motor, low noise, maintenance free, high-speed bearings, smooth operation;
- Two shaker plate sizes, corresponding to different carrying capacity requirements, paired with anti slip silicone pads to ensure that the container is not displaced during operation;
- Four types of pallet clamp, corresponding to different container and sample requirements ;
- LCD display, functional zoning display, intuitive and clear, easy to operate;

Product Model	SK-O270T		
Image			
Power Supply	AC100-240V 50/60Hz		
Power (W)	30		
Display Mode	LCD		
Motion Type	Orbital shaker		
Amplitude (mm)	20		
Speed Range (rpm)	40-300		
Timing Range	1min-99h59min		
Max Load Capacity (with platform)	2kg		
Platform Type	Plate platform	Universal spring rack	Platform with adjust rollers
	 (Only suit 270-type shaker)		