

Specifications

Model	WH260-H
Order No.	W3012603
Maximum temperature setting [°C]	450
Maximum temperature setting with E-sensor [°C]	300
Safety temperature [°C]	50~500 adjustable
High temperature protection [°C]	10~50 adjustable
Temperature stability with E-sensor [°C]	±1
Heat output [W]	800
External temperature sensor	Pt100
PID parameter	3 sets
Speed setting range [rpm]	100~1500
Stirring quantity max [L / H ₂ O]	20
Temperature and Speed display	LCD display
Temperature and Speed setting	Turning knob
Time setting	1~1999min / continuous
IP code	IP21
Top plate material	Ceramic glass
Top plate dimensions [mm]	180×180
RS232 interface	Yes
Dimensions W x D x H [mm]	190x320x85
Weight [kg]	3.1
Power supply	220V/50Hz

*Measured with 500ml water at 80 °C

Standard package includes: Hot plate / stirrer, temperature sensor PT100-01 with holder and clamp, stirrer bar (Order No. 1.230.8)

WIGGENS GmbH

Wiescher Str. 11-13 // Haus 3

42277 Wuppertal Germany

Tel: 0049 202 373 29 58 – 0

WIGGENS China

Room 426, Hall A, Office Building M8, No.1 Jiuxianqiao

East Road, Chaoyang District, Beijing 100015, China

Tel: +86 400-809-2068

Fax: +86 400-809-2068-112

Email: info@wiggins.com service@wiggins.com