

Autoclave Product Description

MODEL	LX-B35L (Digital& hand shaking)	LX-B50L (Digital& hand shaking)	LX-B75L (Digital& hand shaking)	LX-B100L (Digital& hand shaking)
Volume	35L	50L	75L	100L
Max Working Temperature	134°C	134°C	134°C	134°C
Max Working Pressure	0.25Mpa	0.25Mpa	0.25Mpa	0.25Mpa
Time Range	0-9999min	0-9999min	0-9999min	0-9999min
Temperature Range	RT-134°C	RT-134°C	RT-134°C	RT-134°C
Sterilization Chamber	Φ280*477	Φ300*537	Φ360*477	Φ360*622mm
Power	AC220V/50Hz /2.5KW	AC220V/50Hz /2.5KW	AC220V/50Hz /3.5KW	AC220V/50Hz /4.0KW
Drying Temperature	60-110°C	60-110°C	60-110°C	60-110°C

