

Vertical Pressure Steam Sterilizer

Used to medical and health undertakings, scientific research, agriculture and other units, medical equipment, dressing glassware, solution medium disinfection and sterilization.

Production Introduces

- Fully stainless steel structure.
- Hand wheel type of quick-open door structure.
- Door safety lock system.
- Digital display of working status, touch of key.
- Self-inflating type seal.
- Over temperature&over pressure auto-protection.
- Safe protection of water lacking.
- Automatically shut off with beep reminding after sterilization.
- Auto discharge the cool air, and steam discharging automatically after sterilization.
- Used for sterilization of medical instruments, medical cotton products.
- Fully stainless steel SUS304/AISI|304-3MM.
- Possible to install drying system according request.



Technical Parameter

Model	FSF-35HD	FSF-50HD	FSF-75HD	FSF-100HD
Chamber Volume	35L Φ318*415mm	50L Φ340*550mm	75L Φ400*600mm	100L Φ440*650mm
Design Pressure	0.3Mpa			
Design Temperature	150℃			
Max.Working Temperature	138℃			
Max.Working Pressure	0.25MPa			
Timer Range	0-99min or 0-99h 59 min			
Power	1*2.5KW/AC 220V.50HZ	2*1.5KW/AC 220V.50HZ	3*1.5KW/AC 220V.50HZ	
Overall Dimension(mm)	470*450*1020	520*490*1150	550*530*1180	580*560*1290
Transport Dimension(mm)	540*540*1160	580*580*1280	610*610*1320	640*640*1400
G.W/N.W(kg)	84/65	97/78	115/97	135/109