

ACE

GmbH



• Max: 0 ~ 220g
d: 0.0001g



ACE

GmbH
Analytical Balance

Feature:

- Tare function
- Auto zero tracking
- Piececounting function
- Selectable measure units: g, mg, t, ct
- Hook for underbalance weighing
- Selectable filter & speed lever
- 5 inch touch screen
- Internal calibration

A224ABT

Analytical Electronic Balance

Feature :

- Tare function
- Auto zero tracking
- Piececounting function
- Selectable measure units: g, mg, t, ct
- Hook for underbalance weighing
- Selectable filter & speed lever
- 5 inch touch screen display with white backlight
- Internal calibration



Model	A224ABT	
Capacity(g)	220g	
Pan size	Φ80mm	
Sensor	Magnet	
Resolution(g)	0.0001	
Min weighing(g)	0.0004	
Stable Time	≤3s	
Display	5 Inch touch screen	
Material	Aluminum	
Operation temp	5-25℃	
Repeat ability	±0.0002g	
Liner	±0.0003g	
Draft shield size	240mmx190mmx265mm	
Cal.weight	200g	
Option	RS232C & RS485, Printer	
Packing size	Inside: 415x285x430mm	Outside: 500x380x530cm
G . weight	10kg	

ACE

GmbH
Langenhorner Chaussee 600,
22419 Hamburg, Germany

Customer Name :
Company Name :
Address :
Phone :

