

YP Series Electronic Precision Balance (0.01g/0.1g)



YP series **0.1g** specification(LED display, "D" for LCD display)

Model No.	YP2002	YP3002	YP5002	YP6002	YP10002	YP20002	YP30002	YP50002
	YP2002D	YP3002D	YP5002D	YP6002D	YP10002D	YP20002D	YP30002D	YP50002D
	YP2002S	YP3002S	YP5002S	YP6002S	YP10002S	YP20002S	YP30002S	YP50002S
Capacity	0-200g	0-300g	0-500g	0-600g	0-1000g	0-2000g	0-3000g	0-5000g
Readability	0.01g							
Repeatability(≤)	±0.02g							
Linearity(≤)	±0.02g							
Pan Size	125x160mm				143x192mm			
Product Size	22.5x22x8cm							
Dual display Size	27.5x22x8cm							
Packing Size	30.5x24.5x11cm							
N.W.	1200g							
G.W.	1500g							
Power	110V/220-240V AC							



YP series **0.01g** specification (LED display, "D" for LCD display, "S" for dual display)

Model No.	YP2001	YP6001	YP10001	YP15001	YP20001	YP30001	YP50001	YP60001	YP100001
	YP2001D	YP6001D	YP10001D	YP15001D	YP20001D	YP30001D	YP50001D	YP60001D	YP100001D
	YP2001S	YP6001S	YP10001S	YP15001S	YP20001S	YP30001S	YP50001S	YP60001S	YP100001S
Capacity	0-200g	0-600g	0-1000g	0-1500g	0-2000g	0-3000g	0-5000g	0-6000g	0-10000g
Readability	0.1g								
Repeatability(≤)	±0.2g								
Linearity(≤)	±0.2g								
Pan Size	143x192mm								
Product Size	22.5x22x8cm								
Dual display Size	27.5x22x8cm								
Packing Size	30.5x24.5x11cm								
N.W.	1200g								
G.W.	1500g								
Power	110V/220-240V AC								

