



Series Precision Electronic Balance

(1mg/0.01g/0.1g)



High precision strain gauge loadcell, skip the tedious operation process. Dual interface , dual display optional, provide more intelligent applications in the field of strain gauges, and provide customers with good cost performance and convenient use experience.

Main features

- High contrast LCD display
- Big size stainless steel scale pan
- Glass windshield
- Automatic fault detection, overload protection
- Unit conversion, counting, percentage, etc.
- Tare function
- RS232C interface (Option USB interface)
- Below weighing facility (Optional)
- Dual display (Optional)
- 4 x AA battery (Optional)



TF series 1mg specification("S" for dual display model)

Model No.	TF1003	TF2003	TF3003	TF5003
	TF1003S	TF2003S	TF3003S	TF5003S
Capacity	0-100g	0-200g	0-300g	0-500g
Readability	1mg			
Repeatability(≤)	±2mg			
Linearity(≤)	±2mg			
Interface	RS232C/USB(Optional)			
Pan Size	Φ90mm			
Product Size	32x19.5x28cm			
Packing Size	37.5x25.5x34cm			
N.W.	2600g			
G.W.	3500g			
Power	110V/220-240V AC			

