

Transparent Water Bath

◎ Used for precision constant temperature and auxiliary.



Characteristics:

Patent No. : ZL201920769902X

- Using innovative design structure, effectively prevent steam leakage.
- Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.
- Shell, interior and cover of the product are all made of transparent engineering plastics, can clearly see the object to be heated.
- Power off when water shortage to ensure the safe operation of the product.

Main Technical Parameters:

Model	DK-98-IV				
Capacity	9L	2.6L	6L	12L	20L
Voltage (V)	200-240V/100-120V				
Power (kW)	0.8	0.4	0.5	1.0	1.5
Temp. Motion (°C)	±0.5	±0.5	±0.5	±0.5	±0.5
Temp. Range (°C)	Rt+5~100				
Temp. Sensitivity (°C)	≤±1	≤±1	≤±1	≤±1	≤±1
Shelves Height (mm) Distance from bottom	100mm	100mm	100mm	150mm	100mm
Shelves Holes Diameter	Φ 18.7 *48mm	Φ 18.7 *16mm	Φ 18.7 *32mm	Φ 18.7 *48mm	Φ 18.7 *96mm
Inner Chamber Size (L*W*H)	280 x 220 x150	143*130*150	290*145*150	285*225*200	485*285*150
Exterior Size (L*W*H)	400 x260x 190	254*190*193	401*190*193	401*276*242	610*338*200
Packing Size (L*W*H)	450 x350 x270	295*270*275	457*270*275	450*350*300	640*410*230
N.W. (KG)	3.5	2.1	3.1	3.9	6
G.W. (KG)	5	2.7	3.8	4.8	8.5