

Customer Name :
 Company Name :
 Address :
 Phone :

ACE
 GmbH



ACE
 GmbH
 Moisture Analyzer

Feature:
 5inch touch screen
 Real time testing curve
 Standard/Fast/Soft Mode
 Automatic/Manual/Timed stop mode
 Silent fast cooling fan
 Adaptive input voltage
 Test frequency selection
 Clock function
 Historical test report



A100ME
 Economy Moisture analyzer

ACE
 GmbH
 Langenhorner Chaussee 600,
 22419 Hamburg, Germany



Model	A100ME
Capacity	110 g
Readability	0.001 g
Sensor	Load Cell
Pan Size	Φ96mm (20pcs)
Display	5inch touch screen
Operate.Temp	5°C-35°C
Heating Source	Halogen lamp
Temp.Sensor	PT-100
Temp.Range	40°C-160°C
Temp.Setting	1°C
Moisture.Range	100.00%-0.00%
Moisture Readability	0.01%
Testing Range	100.00%-0.00%
DRY% Readability	0.01%
Storage	20
Analysis Curve	Real time dynamic curve
Testing Height	30mm
LANGUAGE	EN
Power supplier	220V 50Hz
Dimension	336X214X166mm
Packing Size	415*310*255mm
G.Weight/N.Weight	4.0kg/3.5kg