

T-MS MULTI SHAKER (products_show-709-174-T_MS_MULTI_SHAKER-INSTRUMENT.html)

|

Topscien Mini Block (products_show-783-174-Topscien_Mini_Block-INSTRUMENT.html)

|

Topscien Incu Block (products_show-782-174-Topscien_Incu_Block-INSTRUMENT.html)

|

Topscien heated lid Incu Block Incubator (products_show-781-174-Topscien_heated_lid_Incu_Block_Incubator-INSTRUMENT.html)

|

Topscien Multi-IncuPlus (products_show-780-174-Topscien_Multi_IncuPlus-INSTRUMENT.html)

|

PGspin(PRF CGF) centrifuge (products_show-745-174-PGspin(PRF_CGF)_centrifuge-INSTRUMENT.html)

|

Digital Water Baths (products_show-779-174-Digital_Water_Baths-INSTRUMENT.html)

|

AccuMini Mini Gradient Thermal Cyclor (products_show-778-174-AccuMini_Mini_Gradient_Thermal_Cycler-INSTRUMENT.html)

|

AccuNeo Gradient Thermal Cyclor (products_show-777-174-AccuNeo_Gradient_Thermal_Cycler-INSTRUMENT.html)

|

THM-INCUB HEATING MICROPLATE INCUBATOR (products_show-708-174-THM_INCUB_HEATING_MICROPLATE_INCUBATOR-INSTRUMENT.html)

|

AccuNeo Plus Gradient Thermal Cyclor (products_show-776-174-AccuNeo_Plus_Gradient_Thermal_Cycler-INSTRUMENT.html)

|

TCM-CAP1A Automatic Capping Machine (products_show-775-174-TCM_CAP1A_Automatic_Capping_Machine-INSTRUMENT.html)

|

Products



([JavaScript:FrontProducts_detail02_scroll\('FrontProducts_detail02-1440567066784_key10ScrollDiv',84\)\)](#))



([JavaScript:FrontProducts_detail02_scroll\('FrontProducts_detail02-1440567066784_key10ScrollDiv',84\)\)](#))

Product details

Spinplus Centrifuge

Spinplus is the most powerful and competitive centrifuge in its class it can be placed 6 of 15ml conical tubes or glass tubes, 8 of 15ml conical tubes or glass tubes conical tubes which with screw caps 1 can be 100- 5000RPM. the LCD display with ability to change between RPM and RCF

Product Name: **Spinplus Centrifuge**

Release time: **2024-10-8 22:07:46**

Company Feature:

TOPSCIEN is a leading provider of OEM or Original Equipment Manufacture liquid handling solutions. Thanks to our know-how and experience, we can take care of the whole process from engineering concept to dispatch. We have developed different technologies that can be applied to custom pipetting solutions. These include both electronic and mechanical pipette platforms. As a customer can rest assured: Topscien's solutions allow you to improve not only the safety and reliability but also the ease-of-use of your liquid handling system.

In many applications, such as diagnostic test and analysis systems, the accuracy and reliability of the system are mission-critical. This requires a lot from all components of the system, including the pipette. Topscien places great emphasis on the accuracy, precision, safety and friendliness of its products. Its application-specific solutions have been carefully manufactured in accordance with international quality standards.

Custom-made for your needs

The Start electronic and mechanical pipettes, as well as the new dispenser module families adapt platforms for your liquid handling system. For more specific needs, Topscien offers tailored solutions based on any of its pipettes.

To learn more about Topscien's OEM solutions, contact us. wilsonw@topscien.com (mailto:wilsonw@topscien.com)

Feature

- Safety lid-lock
- LCD display showing
- Open lid automatically
- The operation interface simple and easy
- 8 or 6 of 15ml conical tubes with screw caps can be placed
- Gentle braking at low speeds with efficient separation acceleration and braking time(20-80 seconds) can be adjustable for the 0-5000rpm speed
- The LCD display with ability to change between RPM andRCF
- Separate short spin button for convenient and quick spins .
- Automatic release the lid when operation has stopped to save processing time
- Operation can be timed from 10 seconds to 99 minutes with bell signal after finish operation
- Safety lid-lock prevent opening the centrifuge lid until the rotor has stopped spinning
- 1 YEAR WARRANTY





Technical Data

Model	TC-Spinplus
Power supply	110v-60Hz or 230v-50/60Hz
Power requirement	40W
Max.rotational speed(rpm)	100-5000rpm (in 100 rpm increments)
Relative centrifugal force(rcf)	3074xg
Timer	15sec-99min
Acceleration time to max.speed	20-80 seconds available
Braking time from max.speed	20-80 seconds available